
Product Data Sheet

Product Name: Belantamab

Cat. No.: GC68369

Chemical Properties

Cas. No. 2061894-48-0

Formula M.Wt

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Belantamab (GSK2857914) is a humanised IgG1 anti-**BCMA (TNFRSF17)** monoclonal antibody. Belantamab can be used in the synthesis of antibody-drug conjugate (ADC), Belantamab mafodotin^[1].

[1]. Ivo Demel, et al. Focus on monoclonal antibodies targeting B-cell maturation antigen (BCMA) in multiple myeloma: update 2021. Br J Haematol. 2021 May;193(4):705-722.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA