
Product Data Sheet

Product Name: Mutant p53 modulator-1

Cat. No.: GC68371

Chemical Properties

Cas. No. 2746371-35-5

Formula $C_{27}H_{32}F_4N_8O_2$ M.Wt 576.59

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Mutant p53 modulator-1 is a mutant p53 modulator. Mutant p53 modulator-1 reduces the progression of cancers that contain a p53 mutation (extracted from patent WO2021231474A1, compound 231B)^[1].

[1]. Binh Vu, et al. METHODS AND COMPOUNDS FOR RESTORING MUTANT p53 FUNCTION. Patent WO2021231474A1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
