

---

**Product Data Sheet**

---

Product Name: FGN849  
Cat. No.: GC68387

**Chemical Properties**

Cas. No. 2069250-01-5

Formula  $C_{42}H_{37}N_5O_6$  M.Wt 707.77

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

**Background**

FGN849 is the active catabolite of IMGN632<sup>[1]</sup>.

[1]. Kovtun Y, et al. A CD123-targeting antibody-drug conjugate, IMGN632, designed to eradicate AML while sparing normal bone marrow cells. Blood Adv. 2018;2(8):848-858.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA