
Product Data Sheet

Product Name: SJ995973

Cat. No.: GC68389

Chemical Properties

Cas. No.

Formula C₃₇H₄₀ClN₇O₅S

M.Wt 730.28

Solubility DMSO : 100 mg/mL (136.93 mM; Need ultrasonic)

Storage 4°C, away from moisture and light

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

SJ995973 (PROTAC) is a uniquely potent degrader of **bromodomain and extra-terminal (BET)** proteins.

[1]. Min J, et al. Phenyl-Glutarimides: Alternative Cereblon Binders for the Design of PROTACs. *Angew Chem Int Ed Engl.* 2021 Dec 13;60(51):26663-26670.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA