
Product Data Sheet

Product Name: 4-Propionamidophenol

Cat. No.: GC68555

Chemical Properties

Cas. No. 1693-37-4

Formula $C_9H_{11}NO_2$

M.Wt 165.19

Solubility DMSO : 250 mg/mL (1513.41 mM; Need ultrasonic)

Storage 4°C, stored under nitrogen

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

4-Propionamidophenol is a p-acetamidophenol analog^[1].

[1]. Chavan, B. B., et al. Identification and characterization of vilazodone metabolites in rats and microsomes by ultrahigh-performance liquid chromatography/quadrupole time-of-flight tandem mass spectrometry. Rapid Communications in Mass Spectrometry. 2017.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
