

Product Data Sheet

Product Name: Barzolvolimab
 Cat. No.: GC68727

Chemical Properties

Cas. No. 2438203-51-9

Formula M.Wt

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Barzolvolimab (CDX 0159) is a humanized anti-**KIT** IgG1 monoclonal antibody. Barzolvolimab specificity and potently inhibits KIT activation by SCF. Barzolvolimab can reduce skin mast cells and disease activity in chronic inducible urticaria^[1].

Barzolvolimab 抑制 KIT 活化 (SCF) ^[1]

[1]. Terhorst-Molawi D, et, al. Anti-KIT antibody, barzolvolimab, reduces skin mast cells and disease activity in chronic inducible urticaria. Allergy. 2022 Nov 16.