

---

## Product Data Sheet

---

Product Name: BCN-PEG4-NHS ester

Cat. No.: GC68745

### Chemical Properties

Cas. No.

Formula C<sub>26</sub>H<sub>38</sub>N<sub>2</sub>O<sub>10</sub>

M.Wt

538.59

Solubility

Storage

Store at -20°C

General tips

For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

### Background

BCN-PEG4-NHS ester is a PEG-based PROTACs<sup>[1]</sup>.

PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins<sup>[1]</sup>.

[1]. "An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562."

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA