
Product Data Sheet

Product Name: Galegenimab

Cat. No.: GC69145

Chemical Properties

Cas. No. 2403683-24-7

Formula M.Wt

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Galegenimab (FHTR 2163; RG 6147; RO 7171009), an anti-**High-temperature requirement A1 (HTRA1)** antibody fragment. Galegenimab can be used for age-related macular degeneration (AMD) research^[1].

[1]. Jonathan B Lin, et al. Neuroprotection for Age-Related Macular Degeneration. *Ophthalmol Sci.* 2022 Jul 5;2(4):100192.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA