

---

## Product Data Sheet

---

Product Name: Ganitumab

Cat. No.: GC69149

### Chemical Properties

Cas. No. 905703-97-1

Formula M.Wt 145.70 kDa

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

Ganitumab (AMG 479) is a recombinant human monoclonal antibody to the **human type 1**

[1]. Moody G, et, al. IGF1R blockade with ganitumab results in systemic effects on the GH-IGF axis in mice. J Endocrinol. 2014 Mar 17;221(1):145-55.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA