
Product Data Sheet

Product Name: GSK3004774

Cat. No.: GC69193

Chemical Properties

Cas. No. 2138814-32-9

Formula $C_{32}H_{40}N_4O_3$ M.Wt 528.69

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

GSK3004774 is a potent, nonabsorbable agonist of **CaSR^[1]**.

[1]. Sparks SM, et al. Identification of potent, nonabsorbable agonists of the calcium-sensing receptor for GI-specific administration. Bioorg Med Chem Lett. 2017 Oct 15;27(20):4673-4677.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
