
Product Data Sheet

Product Name: HO-PEG-CH₂COOH (MW 35000)

Cat. No.: GC69237

Chemical Properties

Cas. No.

Formula

M.Wt 35000

Solubility DMSO : 12.5 mg/mL (0.36 mM; ultrasonic and warming and heat to 60°C)

Store at -20°C,
Storage stored under
nitrogenGeneral
tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

HO-PEG-CH₂COOH (MW 35000) is a **PROTAC molecule**^[1].

[1]. An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
