
Product Data Sheet

Product Name: Human PD-L1 inhibitor I

Cat. No.: GC69243

Chemical Properties

Cas. No. 2135542-86-6

Formula $C_{110}H_{152}N_{26}O_{32}$

M.Wt

2350.58

Solubility

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Human PD-L1 inhibitor I may disturb the binding of hPD-L1 to hPD-1^[1].

[1]. Tianyu Wang, et al. Development of Inhibitors of the Programmed Cell Death-1/Programmed Cell Death-Ligand 1 Signaling Pathway. J Med Chem. 2019 Feb 28;62(4):1715-1730.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA