
Product Data Sheet

Product Name: Lacutamab

Cat. No.: GC69351

Chemical Properties

Cas. No. 2187368-16-5

Formula M.Wt

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Lacutamab (IPH4102) is a humanized monoclonal antibody (mAb) directed against the immune receptor molecule **KIR3DL2**, produced by recombinant technology in CHO cells. Lacutamab can be used in research of cutaneous T-cell lymphoma^[1].

[1]. Van Der Weyden C, et, al. IPH4102, a monoclonal antibody directed against the immune receptor molecule KIR3DL2, for the treatment of cutaneous T-cell lymphoma. Expert Opin Investig Drugs. 2018 Aug;27(8):691-697.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA