

---

## Product Data Sheet

---

Product Name: Lirentelimab

Cat. No.: GC69380

### Chemical Properties

Cas. No. 2283348-97-8

Formula M.Wt

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

Lirentelimab (AK002) is a humanized Completely prevented passive systemic anaplasia (PSA) as shown by a lack of change in rectal temperature and symptom scores in mice.

[1]. Youngblood BA, et al. AK002, a Humanized Sialic Acid-Binding Immunoglobulin-Like Lectin-8 Antibody that Induces Antibody-Dependent Cell-Mediated Cytotoxicity against Human Eosinophils and Inhibits Mast Cell-Mediated Anaplasia in Mice. *Int Arch Allergy Immunol.* 2019;180(2):91-102.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA