

---

## Product Data Sheet

---

Product Name: Lokivetmab

Cat. No.: GC69391

### Chemical Properties

Cas. No. 1533403-95-0

Formula M.Wt

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

Lokivetmab (Anti-Canine IL31 Recombinant Antibody) is an anti-canine **IL-31** monoclonal antibody that can be used for the research of atopic dermatitis (AD) in dogs<sup>[1]</sup>.

[1]. Meichner K, et al. Lokivetmab therapy for pruritus in a dog with cutaneous mastocytosis. Vet Dermatol. 2019 Feb;30(1):73-e22.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA