
Product Data Sheet

Product Name: MLS001006105

Cat. No.: GC69477

Chemical Properties

Cas. No. 455310-36-8

Formula C₁₅H₁₅N₅O₃S₂

M.Wt 377.44

Solubility DMSO : 125 mg/mL (331.18 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

MLS001006105 is used for biological test against the protein target: lamin isoform A-delta10^[1].

[1]. N-[4-(3,4-dimethoxyphenyl)-1,3-thiazol-2-yl]-2-(1H-1,2,4-triazol-5-ylsulfanyl)acetamide

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
