
Product Data Sheet

Product Name: MS9427
Cat. No.: GC69505

Chemical Properties

Cas. No. 2772613-37-1

Formula C₄₈H₅₈ClFN₈O₁₂ M.Wt 993.47

Solubility DMSO : 70 mg/mL (70.46 mM;
Need ultrasonic) Storage 4°C, protect from light, stored
under nitrogen

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it
in the ultrasonic bath for a while. Stock solution can be stored below -20°C for
several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship
with RT , or blue ice upon request.

Structure

Background

MS9427 is a potent **PROTAC**

[1]. Yu X, et al. Exploring Degradation of Mutant and Wild-Type Epidermal Growth Factor Receptors Induced by Proteolysis-Targeting Chimeras. J Med Chem. 2022 Jun 23;65(12):8416-8443.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
