
Product Data Sheet

Product Name: PI3K-IN-30

Cat. No.: GC69707

Chemical Properties

Cas. No. 2281803-22-1

Formula C₂₀H₂₅F₂N₇O₃

M.Wt

449.45

Solubility

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

PI3K-IN-30 (compound 6d) is a potent **PI3Kδ**, respectively^[1].

[1]. Van Dort ME, et al. Structural effects of morpholine replacement in ZSTK474 on Class I PI3K isoform inhibition: Development of novel MEK/PI3K bifunctional inhibitors. Eur J Med Chem. 2022;229:113996.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
