
Product Data Sheet

Product Name: Pomalidomide-5'-PEG8-C2-COOH

Cat. No.: GC69726

Chemical Properties

Cas. No. 2138440-78-3

Formula $C_{32}H_{47}N_3O_{14}$ M.Wt 697.73

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**Pomalidomide-5'-PEG8-C2-COOH is an **PROTACs**^[1].

[1]. Posternak G, et al. Functional characterization of a PROTAC directed against BRAF mutant V600E. Nat Chem Biol. 2020 Nov;16(11):1170-1178.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
