
Product Data Sheet

Product Name: PROTAC CDK12/13 Degradar-1 TFA

Cat. No.: GC69748

Chemical Properties

Cas. No.

Formula C₄₇H₄₇F₃N₁₀O₈

M.Wt 936.93

Solubility DMSO : 25 mg/mL (26.68 mM; ultrasonic and warming and heat to 60°C)

Storage Store at -20°C, away from moisture and light

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

-

[1]. Jianzhang Yang, et al. Discovery of a Highly Potent and Selective Dual PROTAC Degradar of CDK12 and CDK13. J Med Chem. 2022 Aug 8.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA