
Product Data Sheet

Product Name: Quetmolimab

Cat. No.: GC69789

Chemical Properties

Cas. No. 2084037-83-0

Formula M.Wt 150 kDa

Solubility Storage Store at -80°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Quetmolimab is a humanized anti-**CX3CL1** is a chemokine with a modulating effect on chemotaxis and adhesion^[1].

[1]. Blaess J, et al. Immunosuppressive agents for rheumatoid arthritis: a systematic review of clinical trials and their current development stage. Ther Adv Musculoskelet Dis. 2020 Dec 16;12:1759720X20959971.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA