
Product Data Sheet

Product Name: Rocatinlimab

Cat. No.: GC69826

Chemical Properties

Cas. No. 2431972-52-8

Formula M.Wt

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Rocatinlimab (AMG 451) (KHK4083) is a fully human, non-fucosylated, **immunoglobulin G1 (IgG1) anti-OX40 monoclonal antibody**. Rocatinlimab can be used for the research of atopic dermatitis (AD)^[1].

[1]. Ana Maria LÉ, et al. OX40-OX40L Inhibition for the Treatment of Atopic Dermatitis-Focus on Rocatinlimab and Amlitelimab. *Pharmaceutics*. 2022 Dec 8;14(12):2753.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA