
Product Data Sheet

Product Name: Sibeprenlimab

Cat. No.: GC69887

Chemical Properties

Cas. No. 2382896-07-1

Formula M.Wt

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Sibeprenlimab (VIS649) is a **humanized IgA nephropat (IgAN)[1]**.

[1]. Mohit Mathur, et al. Safety, Tolerability, Pharmacokinetics, and Pharmacodynamics of VIS649 (Sibeprenlimab), an APRIL-Neutralizing IgG2 Monoclonal Antibody, in Healt Volunteers. Kidney Int Rep. 2022 Feb 8;7(5):993-1003.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA