
Product Data Sheet

Product Name: Tomaralimab

Cat. No.: GC70042

Chemical Properties

Cas. No. 1449294-76-1

Formula M.Wt 145.32 kDa

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Tomaralimab (OPN-305) is a humanised **anti-TLR2** IgG4 monoclonal antibody.

Tomaralimab has the potential for the research of Myelodysplastic Syndromes (MDS)^[1].

[1]. Guillermo Garcia-Manero MD, et al. A Clinical Study of Tomaralimab (OPN-305), a Toll-like Receptor 2 (TLR-2) Antibody, in Heavily Pre-Treated Transfusion Dependent Patients with Lower Risk Myelodysplastic Syndromes (MDS) That Have Received and Failed on Prior pometlating Agent (HMA) Therapy. Blood, 2018, 798.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA