
Product Data Sheet

Product Name: Tozorakimab

Cat. No.: GC70045

Chemical Properties

Cas. No. 2376858-66-9

Formula CH₃ M.Wt 15.03

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Tozorakimab (MEDI-3506) is a human immunoglobulin G1 monoclonal antibody targeting **interleukin-33**. Tozorakimab can be used to research chronic obstructive pulmonary disease^[1].

[1]. F. Reid, et al. Safety, Pharmacokinetics and Immunogenicity of the Anti-Interleukin-33 Antibody, Tozorakimab: Results of a Phase 1 Study in Healthy Adults and Patients with Chronic Obstructive Pulmonary Disease.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA