
Product Data Sheet

Product Name: Zagotenemab

Cat. No.: GC70163

Chemical Properties

Cas. No. 2019133-28-7

Formula M.Wt

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Zagotenemab (LY3303560) is a humanised anti-**tau** antibody that selectively binds and neutralises **tau** deposits in the brain. Zagotenemab can be used in Alzheimer's disease research^[1].

[1]. Abyadeh M, et al. Comparative Analysis of Aducanumab, Zagotenemab and Pioglitazone as Targeted Treatment Strategies for Alzheimer's Disease. Aging Dis. 2021 Dec 1;12(8):1964-1976.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA