
Product Data Sheet

Product Name: Anti-Mouse CD106/VCAM-1 Antibody (M/K-2.7)

Cat. No.: GC72330

Chemical Properties

Cas. No.

Formula M.Wt 150000

Solubility Storage Store at 4'C, do not freeze.

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Anti-Mouse CD106/VCAM-1 Antibody (M/K-2.7) is a rat-derived IgG1, κ type antibody inhibitor, targeting to mouse CD106/VCAM-1.

The antibody framework is stable, specific and adaptable, and has the ability to bind both antigens and endogenous immune receptors. Monoclonal antibodies have several derivatives, including bispecific antibodies, antibody-drug conjugates, and antibody fragments, and have significant effects in fields such as immunology and oncology. When designing inhibitory antibodies, considerations include identification of antigen-specific variable regions, choice of expression system, use of multispecific formats, and antibody derivatives based on fragmentation, oligomerization, or conjugation with other functional moieties^[1].

<p >References:

[1]. Goulet DR, Atkins WM. Considerations for the Design of Antibody-Based Therapeutics. J Pharm Sci. 2020 Jan;109(1):74-103.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA