
Product Data Sheet

Product Name: Anti-Mouse 4-1BB/CD137 Antibody (3H3)

Cat. No.: GC72337

Chemical Properties

Cas. No.

Formula M.Wt 150000

Solubility Storage 4°C, do not freeze.

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Anti-Mouse 4-1BB/CD137 Antibody (3H3) is an anti-mouse 4-1BB IgG2a antibody inhibitor derived from the host Rat.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA