
Product Data Sheet

Product Name: Cannabidiolic Acid-d3 (CRM)
Cat. No.: GC91393

Chemical Properties

Cas. No.

Formula $C_{22}H_{27}D_3O_4$ M.Wt 361.5

Solubility Storage Store at $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Cannabidiolic acid-d3 GC LC-MS