

Product Data Sheet

Product Name: PAMAM Dendrimer G2.0 Amine (water solution)

Cat. No.: GC91476

Chemical Properties

Cas. No. 26937-01-9

Formula NH2(CH2)2NH2]:(G=2);dendri PAMAM(NH2)16 M.Wt 3256.2

Solubility Water: Soluble Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

PAMAM G2.0 PAMAM 29 16 PAMAM G2.0 15 / 100 500 HEK293T HeLa 700 PAMAM G2.0

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA