

Product Data Sheet

Product Name: 2×Taq PCR Super Fast Mix(with blue dye)
Cat. No.: GK10035

Features

Applications	1. 基因组DNA、cDNA、质粒DNA、PCR产物、DNA模板、2. 检测PCR产物、SNP、3. PCR产物3'端dA、TA
Shipping	常温
Storage	-20°C保存24个月、4°C保存6个月
Usage	For Research Use Only! Not For Use in Humans.

Protocol

1. 基因组DNA、cDNA、质粒DNA、PCR产物、DNA模板

2×Taq PCR Super Fast Mix(with blue dye)

基因组DNA 10~1000 ng、cDNA 1~30 ng、RT-PCR产物 1~2μl、PCR产物、DNA模板、cDNA、PCR产物、SNP、PCR产物3'端dA、TA

FAQs:

1. 基因组DNA、cDNA、质粒DNA、PCR产物、DNA模板、2. 检测PCR产物、SNP、3. PCR产物3'端dA、TA

PCR产物、SNP、PCR产物3'端dA、TA

基因组DNA、cDNA、质粒DNA、PCR产物、DNA模板、2×Taq PCR Fusion-Fast Mix(with blue dye) (GK10036)

PCR产物、SNP、PCR产物3'端dA、TA

Caution: Product has not been fully validated for medical applications. For research use only.
Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com
Address: 10292 Central Ave. #205, Montclair, CA, USA

