
Product Data Sheet

Product Name: Presapogenin CP4

Cat. No.: GN10329

Chemical Properties

Cas. No. 75799-18-7

Formula C₄₆H₇₄O₁₅

M.Wt

867.07

Solubility

Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Extracted from Anemone hupehensis; Store the product in sealed, cool and dry condition

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA