

---

## Product Data Sheet

---

Product Name: FGF 1 Mouse  
 Cat. No.: GP20206  
 Batch No.: 1

### Product Data

Purity	>98%	Source	Escherichia Coli.
Physical Appearance	solid	Shipping Condition	Shipped at Room temp.
Synonyms	HBGF-1; ECGF-beta; FIBP; FGFBP; FIBP-1; ECGF; ECGFA; GLIO703; FGF1; FGF-a.		
Amino Acid Sequence	MFNLPLGNYK KPALLYCSNG GHFLRILPDG TVDGTRDRSD QHIQLQLSAE SAGEVYIKGT ETGQYLAMDT EGLLYGSQTP NEECLFLERL EENHYNTYTS KKHAEKNWFV GLKKNNGSCKR GPRTHYGQKA ILFLPLPVSS D.		
Solubility	It is recommended to reconstitute the lyophilized Fibroblast Growth Factor-acidic in sterile 18MΩ-cm H <sub>2</sub> O not less than 100μg/ml, which can then be further diluted to other aqueous solutions.		
Formulation	Lyophilized from a 0.2μm filtered concentrated solution in PBS, pH 7.5.		

### Introduction

Acidic fibroblast growth factor is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities, and are involved in a variety of biological processes, including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth and invasion. This protein functions as a modifier of endothelial cell migration and proliferation, as well as an angiogenic factor. It acts as a mitogen for a variety of mesoderm- and neuroectoderm-derived cells in vitro, thus is thought to be involved in organogenesis. Three alternatively spliced variants encoding different isoforms have been described. The binding growth factors are angiogenic agents in vivo and are potent mitogens for a variety of cell types in vitro. There are differences in the tissue distribution and concentration of these 2 growth factors.

### Biological Activity

The ED50 as determined by the dose-dependent proliferation of mouse BALB/c 3T3 cells,

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---

## Product Data Sheet

---

is less than 0.2ng/ml corresponding to a Specific Activity of 5x10<sup>6</sup> IU/mg.

### Stability

Lyophilized Fibroblast Growth Factor-1 although stable at room temperature for 3 weeks, should be stored desiccated below -18°C . Upon reconstitution FGF-a should be stored at 4°C between 2-7 days and for future use below -18°C .Please prevent freeze-thaw cycles.

### Background

Fibroblast Growth Factor-acidic Mouse Recombinant (FGF-1) produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 141 amino acids and having a molecular mass of 15.9kDa.The FGF acidic is purified by proprietary chromatographic techniques.

**Caution: Product has not been fully validated for medical applications. For research use only.**

**Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com**

**Address: 10292 Central Ave. #205, Montclair, CA, USA**