
Product Data Sheet

Product Name: IL 15 Mouse, His
 Cat. No.: GP20572
 Batch No.: 1

Product Data

Purity	>98%	Source	Escherichia Coli.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Synonyms	IL-15; MGC9721; Interleukin-15; IL15; AI503618.		
Amino Acid Sequence	MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSHMNW IDVRYDLEKI ESLIQSIHID TTYTDSDFH PSCKVTAMNC FLLELQVILHEYSNMTLNET VRNVLYLANS TLSSNKNVAE SGCKECEELE EKTFTFLQS FIRIVQMFIN TS.		
Formulation	The mouse IL-15 protein solution contains PBS, pH-7.4 and 10% glycerol.		

Introduction

The protein encoded by this gene is a cytokine that regulates T and natural killer cell activation and proliferation. This cytokine and interleukine 2 share many biological activities. They are found to bind common hematopoietin receptor subunits, and may compete for the same receptor, and thus negatively regulate each other's activity. The number of CD8+ memory cells is shown to be controlled by a balance between this cytokine and IL2. This cytokine induces the activation of JAK kinases, as well as the phosphorylation and activation of transcription activators STAT3, STAT5, and STAT6. Studies of the mouse counterpart suggested that this cytokine may increase the expression of apoptosis inhibitor BCL2L1/BCL-x(L), possibly through the transcription activation activity of STAT6, and thus prevent apoptosis. Two alternatively spliced transcript variants of this gene encoding the same protein have been reported.

Biological Activity

Measured in a cell proliferation assay using CTLL2 mouse cytotoxic T cells. The ED50 for this effect is < 1 ng/ml, corresponding to a specific activity of 1,000,000IU/mg.

Stability

Store at 4°C if entire vial will be used within 2-4 weeks. Store. frozen at -20°C for longer

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Background

Interleukin-15 Mouse Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 152 amino acids (49-162 a.a.) and having a molecular mass of 17.6 kDa. The IL-15 is fused to a 37 amino acid His Tag at N-terminus and purified by proprietary chromatographic techniques.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA