
Product Data Sheet

Product Name: IL 17 Human, His
Cat. No.: GP20581
Batch No.: 1

Product Data

| | | | |
|---------------------|--|--------------------|-------------------------|
| Purity | >98% | Source | Escherichia Coli. |
| Physical Appearance | solid | Shipping Condition | Shipped with Ice Packs. |
| Synonyms | CTLA-8; IL-17; IL-17A; Cytotoxic T-lymphocyte-associated antigen 8; IL17A. | | |
| Formulation | Interleukin-17 protein solution (0.171mg/ml) is supplied in 25mM Na-Acetate pH 4.8 and 50% glycerol. | | |

Introduction

IL17 is a proinflammatory cytokine produced by activated T cells. IL-17 regulates the activities of NF-kappaB and mitogen-activated protein kinases. Interleukin-17 can stimulate the expression of IL6 and cyclooxygenase-2 (PTGS2/COX-2), as well as enhance the production of nitric oxide (NO). High levels of IL-17 are associated with several chronic inflammatory diseases including rheumatoid arthritis, psoriasis and multiple sclerosis.

Stability

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Please avoid freeze thaw cycles.

Background

Interleukin-17A Human Recombinant produced in E.Coli is a single, non-glycosylated, Polypeptide chain containing 132 amino acids fragment (24-155) having a molecular weight of 19.62kDa and fused with a 4.5kDa amino-terminal hexahistidine tag. The IL-17A His is purified by proprietary chromatographic techniques.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA