
Product Data Sheet

Product Name: IL 24 Human
 Cat. No.: GP20629
 Batch No.: 1

Product Data

Purity	>98%	Source	Sacharomyces cerevisiae.
Physical Appearance	solid	Shipping Condition	Shipped at Room temp.
Synonyms	C49A; FISP; MDA7; ST16; IL-24; IL10B; Mob-5; MDA-7; Suppression of tumorigenicity 16 protein; Melanoma differentiation-associated gene 7 protein.		
Solubility	It is recommended to reconstitute the lyophilized IL-24 in sterile 18MΩ-cm H ₂ O not less than 100μg/ml, which can then be further diluted to other aqueous solutions.		
Formulation	Lyophilized from a 0.2μm filtered solution in PBS with BSA as a carrier.		

Introduction

IL24 is a member of the IL10 family of cytokines. It was identified as a gene induced during terminal differentiation in melanoma cells. IL-10B encoded can induce apoptosis selectively in various cancer cells. Overexpression IL-24 leads to elevated expression of several GADD family genes, which correlates with the induction of apoptosis. The phosphorylation of mitogen-activated protein kinase 14 (MAPK7/P38), and heat shock 27kDa protein 1 (HSPB2/HSP27) are found to be induced by this gene in melanoma cells, but not in normal immortal melanocytes. Alternatively spliced transcript variants encoding distinct isoforms have been reported. The glycosylation is essential for activity of IL-24. Functionally, IL-24 has diverse activities. At low concentrations, it induces type I proinflammatory cytokines such as IFN-g, IL-1b, IL-12 and TNF-a. At high concentration, it is a strong inducer of apoptosis in tumor cells, but not normal cells. mda-7/IL-24 is being hailed as a 6omagic bullet6n for cancer gene therapy.

Biological Activity

Measured by its ability to bind to the cell receptor of Capan-1 cells line resulted in Stat-3 activation. The ED50 for this effect is typically 1.0 ng/ml, corresponding to a Specific

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Activity of 1x10⁶ units/mg.

Stability

Lyophilized MDA7 although stable at room temperature for 3 weeks, should be stored desiccated below -18°C . Upon reconstitution MDA-7 Recombinant should be stored at 4°C between 2-7 days and for future use below -18°C .For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA).Please prevent freeze-thaw cycles.

Background

Interleukin 24 human recombinant produced in yeast is a single, glycosylated, polypeptide chain containing 158 amino acids and having a molecular mass of 18 kDa. As a result of glycosylation, the protein migrates at 19.5 kDa on SDS-PAGE.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA