
Product Data Sheet

Product Name: LHRH Human
 Cat. No.: GP21275
 Batch No.: 1

Product Data

Purity	>98%	Source	
Physical Appearance	solid	Shipping Condition	Shipped at Room temp.
Synonyms	Progonadoliberin-1; Progonadoliberin I; LHRH; GRH; GNRH; LNRH.		
Amino Acid Sequence	Pyr-His-Trp-Ser-Tyr-Gly-Leu-Arg-Pro-Gly-NH ₂ .		
Solubility	It is recommended to reconstitute the lyophilized LHRH in sterile 18MΩ-cm H ₂ O not less than 100μg/ml, which can then be further diluted to other aqueous solutions. The LHRH is also soluble in a 1% glacial acetic acid.		
Formulation	The LHRH was lyophilized from a concentrated (1 mg/ml) solution with no additives.		

Introduction

Gonadotropin-releasing hormone 1 (GNRH1), also known as Luteinising-hormone releasing hormone (LHRH), is a peptide hormone responsible for the release of FSH and LH from the anterior pituitary. GNRH1 is synthesized and released by the hypothalamus. At the pituitary, GNRH1 stimulates the synthesis and secretion of the gonadotropins follicle-stimulating hormone (FSH) and luteinizing hormone (LH). These processes are controlled by the size and frequency of GNRH1 pulses, as well as by feedback from androgens and estrogens. Low frequency GNRH1 pulses lead to FSH release, whereas high frequency GNRH1 pulses stimulate LH release. There are differences in GNRH1 secretion between males and females. In males, GNRH1 is secreted in pulses at a constant frequency, but in females the frequency of the pulses varies during the menstrual cycle and there is a large surge of GNRH1 just before ovulation. GNRH1 secretion is pulsatile in all vertebrates, and is necessary for correct reproductive function. Thus, a single hormone, GNRH1, controls a complex process of follicular growth, ovulation, and corpus luteum maintenance in the female, and

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

spermatogenesis in the male.

Stability

Lyophilized LHRH although stable at room temperature for 3 weeks, should be stored desiccated below -18°C . Upon reconstitution LNRH should be stored at 4°C between 2-7 days and for future use below -18°C . For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Background

Lutenizing Hormone Releasing Hormone Human Synthetic is a single, non-glycosylated, polypeptide chain containing 10 amino acids and having a molecular mass of 1182.2 Dalton. The molecular formula is $\text{C}_{55}\text{H}_{75}\text{N}_{17}\text{O}_{13}\cdot\text{C}_2\text{H}_4\text{O}_2$. The CAS Number is 71447-49-9.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA