
Product Data Sheet

Product Name: BCOADC-E2 Human

Cat. No.: GP21436

Batch No.: 1

Product Data

Purity >98%

Source Sf9 insect cells.

Physical
Appearance solidShipping
Condition Shipped with Ice Packs.

Formulation DLST (0.91mg/ml) is supplied in 20mM HEPES buffer pH-8.0, 200mM NaCl and 20% glycerol.

Stability

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.

Background

BCOADC-E2 is a Full-length cDNA coding for the mature form of the human BCOADC-E2 protein having a molecular mass of 47,321 Dalton (pH 6.5). BCOADC-E2 protein is fused to a hexa-histidine purification tag.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA