
Product Data Sheet

Product Name: CTSD Human
Cat. No.: GP21504
Batch No.: 1

Product Data

Purity	>98%	Source	Escherichia Coli.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Synonyms	Cathepsin D; EC 3.4.23.5; CTSD; CPSD; CLN10; MGC2311.		
Formulation	Cathepsin D at 0.1mg/ml in 50mM Tris-Acetate, pH7.5, 1mM EDTA and 20% Glycerol.		

Introduction

Cathepsin D is synthesized as a 54kDa precursor, which is proteolytically processed to an intermediate 48kDa single chain, which matures into more stable 34kDa and 14kDa two chain form. It is an estrogen-regulated lysosomal protease that has been suggested to facilitate cancer cell migration and invasion by digesting the basement membrane, extracellular matrix, and connective tissue. Because of its mitogenic and proteolytic activities, it has been implicated as a prognostic marker in many tumor types. Cathepsin D is expressed in epithelial cells as well as in macrophages.

Stability

Store vial at -20°C to -80°C . When stored at the recommended temperature, this protein is stable for 12 months. Please prevent freeze-thaw cycles.

Background

Cathepsin D Human Recombinant full length protein expressed in E.coli, shows a 48 kDa band on SDS-PAGE. The Cathepsin D is purified by proprietary chromatographic techniques.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA