
Product Data Sheet

Product Name: Enterokinase Porcine

Cat. No.: GP21607

Batch No.: 1

Product Data

Purity	>98%	Source	Porcine.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Synonyms	Enteropeptidase; EC 3.4.21.9; Enterokinase; Serine protease 7; ENTK; MGC133046.		
Formulation	2 IU/μl, 50mM Tris-HCl, pH 8.0, 0.5M NaCl and 50% glycerol.		

Introduction

Enteropeptidase or enterokinase is an enzyme involved in human digestion. It is produced by cells in the duodenum wall, and is secreted from duodenum's glands, the crypts of Lieberk^hn, whenever ingested food enters the duodenum from the stomach. Enteropeptidase has the critical job of turning trypsinogen (a zymogen) to trypsin, indirectly activating a number of pancreatic digestive enzymes. Enteropeptidase is a serine protease enzyme (EC 3.4.21.9). Enteropeptidase is a part of the Chymotrypsin-clan of serine proteases, and is structurally similar to these proteins.

Stability

One year when stored at -20°C, one week at room temperature.

Background

Porcine enteropeptidase is a specific protease which cleaves after the lysine at its recognition site: Asp-Asp-Asp-Asp-Lys. Enterokinase will not cleave a site followed by proline. Theoretical Mw is 21,880 Dalton, the apparent Mw on SDS-PAGE is about 40 kDa. If a fusion tag is located in the N-terminus with an enterokinase site, enterokinase will be able to remove the fusion tag and to generate the protein exactly as you need without adding any unwanted residues. Glpbio's enterokinase is a highly purified enterokinase from porcine. The enzyme has been extensively purified and tested to ensure that there are no other contaminating proteases.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA