

---

**Product Data Sheet**

---

Product Name: AFP Human  
Cat. No.: GP22684  
Batch No.: 1

**Product Data**

Purity	>98%	Source	Human cord serum.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Synonyms	Alpha-fetoprotein; Alpha-fetoglobulin; Alpha-1-fetoprotein; AFP; FETA; HPAFP.		
Formulation	AFP protein filtered (0.2µm) solution in Tris buffered saline pH 7.5 and less than 0.1% NaN <sub>3</sub> .		

**Introduction**

AFP is normally synthesized in the liver, intestinal tract, and yolk sac of the fetus. Antibody to AFP has been shown to be useful in detecting hepatocellular carcinomas (HCC) and germ cell neoplasms, especially yolk sac tumors.

**Stability**

Human Alpha-Fetoprotein should be stored at 2-8°C .

**Background**

Human alpha-fetoprotein purified from pooled human cord serum.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA