
Product Data Sheet

Product Name: ECGS Bovine
Cat. No.: GP22868
Batch No.: 1

Product Data

Purity		Source	Bovine Hypothalamus.
Physical Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.	Shipping Condition	Shipped at Room temp.
Solubility	water: soluble (1 vial reconstituted with 1 ml 0.85% NaCl)		
Formulation	The Endothelial-Mitogen was lyophilized from a concentrated solution containing 0.1M NaCl and 5mM NaH ₂ PO ₄ , pH-7.4.		

Stability

2-8°C

Background

Endothelial Mitogen can be used for Human & Mammalian Vascular Endothelial Cells at concentration of 0.1-0.3mg/ml of serum supplemented media & 0.01-0.1mg with heparin. Endothelial Mitogen can be used for Hybridoma Cell Cloning at concentrations of 0.025-0.15mg/ml of serum supplemented media. Endothelial Mitogen can be used for Balb/C-3T3 Cells at concentrations of 0.2mg/ml of serum Supplemented media.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA