

---

## Product Data Sheet

---

Product Name: BP-1 Human  
Cat. No.: GP22892  
Batch No.: 1

### Product Data

Purity	>98%	Source	Escherichia Coli.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Synonyms	Homeobox protein DLX-4; DLX-7; DLX-8; Beta protein 1; BP1; DLX7; DLX8; DLX9; DLX4.		
Formulation	BP-1 protein in 50mM Tris-Acetate, pH7.5, 1mM EDTA and 20% Glycerol.		

### Introduction

BP-1 gene, a member of the DLX homeobox gene family, is overexpressed in 80% of breast cancer patients, and thus represents a potential new target for diagnosis and treatment. Members of the Dlx gene family contain a homeobox that is related to that of Distal-less (Dll), a gene expressed in the head and limbs of the developing fruit fly (*Drosophila*). The Distal-less (Dlx) family of genes comprises at least 6 different members, DLX1-DLX6. The DLX proteins are postulated to play a role in forebrain and craniofacial development. Three transcript variants have been described for this gene, however, the full length nature of one variant has not been described. Studies of the two splice variants revealed that one encoded isoform functions as a repressor of the beta-globin gene while the other isoform lacks that function.

### Stability

Store vial at -20°C to -80°C . When stored at the recommended temperature, this protein is stable for 12 months. Please prevent freeze-thaw cycles.

### Background

BP-1 Human Recombinant (aa 103-168) expressed in *E. coli*, shows a 35 kDa band on SDS-PAGE. The BP-1 is purified by proprietary chromatographic techniques.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA