
Product Data Sheet

Product Name: BPI Human
Cat. No.: GP22893
Batch No.: 1

Product Data

Purity	>98%	Source	
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Synonyms	Bactericidal/permeability-increasing protein; CAP 57.		
Formulation	BPI (1.15mg/ml) is supplied in 25mM Sodium Acetate pH-4.6.		

Introduction

BPI encodes a lipopolysaccharide binding protein. BPI is related to human neutrophil granules and has bactericidal activity on gram-negative organisms.

Stability

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.

Background

The native protein BPI is isolated from human buffy coats and has a molecular mass of 53kDa.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA