

---

## Product Data Sheet

---

Product Name: C4c Human  
Cat. No.: GP22940  
Batch No.: 1

### Product Data

Purity	>98%	Source	Human Plasma.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Synonyms	Complement C4c; Complement Component C4c; C4c.		
Formulation	The Human Complement C4c solution is in Sodium Phosphate buffer pH 7.0 containing 0.15M NaCl and 0.09% NaN <sub>3</sub> .		

### Introduction

Complement C4c is a degradation product of C4b which is cleaved by C4b/C3b inactivator to yield C4d and C4c. Complement C4c interacts with C1 and C2 to form C3 convertase of the classic activation pathway. Systemic lupus erythematosus (SLE) is often associated with congenital C4 deficiency. Reduced levels of components of the classic pathway (C1, C2, C4, C3) are common after such activation, e.g. in SLE, acute serum sickness and conditions associated with circulating immune complexes.

### Stability

Human C4c although stable at 4°C for 1 week, should be stored at -15°C. Please avoid freeze-thaw cycles.

### Background

Human C4c produced in Human plasma having a molecular mass of 210 kDa.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA