
Product Data Sheet

Product Name: CD19 Human
 Cat. No.: GP23019
 Batch No.: 1

Product Data

Purity	>98%	Source	Sf9, Baculovirus cells.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Synonyms	CD19 Molecule; B-Lymphocyte Surface Antigen B4; T-Cell Surface Antigen Leu-12; Differentiation Antigen CD19; CD19 Antigen; B-Lymphocyte Antigen CD19; CVID3; B4.		
Amino Acid Sequence	EEPLVVKVEE GDNAVLQCLK GTSDGPTQQL TWSRESPLKP FLKLSLGLPG LGIHMRLPLAI WLFIFNVSQQ MGGFYLCQPG PPSEKAWQPG WTVNVEGSGE LFRWNVSDLG GLGCGLKNRS SEGPSSPSGK LMSPKLYVWA KDRPEIWEGE PPCLPPRDSL NQSLSQDLTM APGSTLWLSC GVPPDSVSRG PLSWTHVHPK GPKSLLSLEL KDDRPARDMW VMETGLLLPR ATAQDAGKYY CHRGNLTMSF HLEITARPVL WHWLLRTGGW KLEPKSCDKT HTCPPCPAPE LLGGPSVFLF PPKPKDTLMI SRTPEVTCVV VDVSHEDPEV KFNWYVDGVE VHNAKTKPRE EQYNSTYRVV SVLTVLHQDW LNGKEYKCKV SNKALPAPIE KTISKAKGQP REPQVYTLPP SRDELTKNQV SLTCLVKGFY PSDIAVEWES NGQPENNYKT TPPVLDSGDS FFLYSKLTVD KSRWQQGNVF SCSVMHEALH NHYTQKSLSL SPGKHHHHHH		
Formulation	CD19 protein solution (0.25mg/ml) contains Phosphate Buffered Saline (pH 7.4) and 10% glycerol.		

Introduction

CD19, also known as B-Lymphocyte Surface Antigen B4, belongs to the Ig superfamily which is expressed on the surface of all B-lymphoid cells with the exception of terminally differentiated plasma cells. CD19 is involved as a signal-transducing receptor in the control of differentiation as well as proliferation. CD19 is as an adaptor protein which drafts cytoplasmic signaling proteins to the membrane and operates within the CD19/CD21 complex to decrease the threshold for B cell receptor signaling pathways. CD19 is present on all B cells, and therefore is a biomarker for lymphoma diagnosis, B

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

lymphocyte development and can be used as a target for leukemia immunotherapies.

Stability

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Background

CD19 produced in Sf9 Baculovirus cells is a single, glycosylated polypeptide chain containing 510 amino acids (21-291a.a.) and having a molecular mass of 57.0kDa. (Molecular size on SDS-PAGE will appear at approximately 50-70kDa). CD19 is expressed with a 239 amino acid hIgG-His tag at C-Terminus and purified by proprietary chromatographic techniques.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA