
Product Data Sheet

Product Name: CD34 Human
 Cat. No.: GP23050
 Batch No.: 1

Product Data

Purity	>98%	Source	E.coli.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Synonyms	CD34 molecule; Hematopoietic progenitor cell antigen CD34; CD34 Antigen; CD34.		
Amino Acid Sequence	MGSSHHHHHH SSGLVPRGSH MSLDNNGTAT PELPTQGTFS NVSTNVSQYE TTPSTLGST SLHPVSQHGN EATTNITETT VKFTSTSVIT SVYGNTNSSV QSQTSVISTV FTTPANVSTP ETTLKPSLSP GNVSDLSTTS TSLATSPTKP YTSSSPILSD IKAEIKCSGI REVKLTQGIC LEQNKTSSCA EFKKDRGEGE ARVLCGEEQA DADAGAQVCS LLLAQSEVRP QCLLLVLANR TEISSKLQLM KKHQSDLKKL GILDFTEQDV ASHQSYSQKT.		
Formulation	CD34 protein solution (0.25mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 0.1M NaCl, 20% glycerol and 1mM DTT.		

Introduction

CD34 is a single chain transmembrane glycoprotein which is selectively expressed on human lymphoid and myeloid hemapoietic progenitor cells. CD34 has been used to measure angiogenesis, which reportedly predicts tumor recurrence.

Stability

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Background

CD34 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 280 amino acids (32-290a.a.) and having a molecular mass of 29.7kDa. CD34 is fused to a 21 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA