
Product Data Sheet

Product Name: CD34 Human, GST
Cat. No.: GP23052
Batch No.: 1

Product Data

Purity	>98%	Source	Escherichia Coli.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Synonyms	Hematopoietic progenitor cell antigen CD34; CD34.		
Formulation	CD34 in 50mM Tris-HCl, pH7.5, 10mM L-glutathione (reduced).		

Introduction

CD34 is a single chain transmembrane glycoprotein which is selectively expressed on human lymphoid and myeloid hemapoietic progenitor cells. CD34 has been used to measure angiogenesis, which reportedly predicts tumor recurrence.

Stability

Store vial at -20°C to -80°C . When stored at the recommended temperature, this protein is stable for 12 months. Please prevent freeze-thaw cycles.

Background

CD34 Human Recombinant full length expressed in E.coli, CD34 is antibody reactive. The CD34 is purified by proprietary chromatographic techniques.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA