
Product Data Sheet

Product Name: CST3 Human

Cat. No.: GP23251

Batch No.: 1

Product Data

Purity	>98%	Source	Human urine.
Physical Appearance	solid	Shipping Condition	Shipped at Room temp.
Synonyms	Cystatin-C; Cystatin-3; Neuroendocrine basic polypeptide; Gamma-trace; Post-gamma-globulin; CST3; MGC117328.		
Solubility	Use of a 0.02M NH ₄ HCO ₃ buffer is recommended.		
Formulation	Lyophilized from 0.02M NH ₄ HCO ₃ . May contain traces of buffer salts.		

Introduction

Cystatins are a superfamily of cysteine proteinase inhibitors found in both plants and animals. They comprise a group of proteinase inhibitors, widely distributed in tissues and body fluids, and form tight complexes with cysteine proteases such as cathepsin B, H, L and S. Cystatin C, a secreted molecule of this family, is of interest from biochemical, medicine and evolutionary points of view. Cystatin C, with molecular weight of 13260 Da, is composed of 120 amino acids, lacks carbohydrate and has two disulfide bridges located near the carboxyl terminus. Cystatin C is increased in patients with malignant diseases, and is related to the insufficiency of renal function and appears to be a better marker than creatinine. On the other hand, low levels of cystatin C involve cause the breakdown of the elastic laminae and, subsequently, the atherosclerosis and abdominal aortic aneurysm.

Stability

Human Cystatin-C although stable at 4°C for 1 week, should be stored at -15°C .Please avoid freeze-thaw cycles.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA