
Product Data Sheet

Product Name: MRFAP1L1 Human
 Cat. No.: GP23949
 Batch No.: 1

Product Data

Purity	>98%	Source	Escherichia Coli.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Synonyms	Morf4 family associated protein 1-like 1; PP784; MRFAP1L1.		
Amino Acid Sequence	MGSSHHHHHH SGLVPRGSH MGSMRPLDID EVEAPEEVEV LEPEEDFEQF LLPVINEMRE DIASLIREHG RAYLRTRSKL WEMDNMLIQI KTQVEASEES ALNHVQHPSG EADERVSELC EKAEEKAKEI AKMAEMLVEL VWRIERSESS.		
Formulation	MRFAP1L1 protein solution (0.5mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 0.15M NaCl, 20% glycerol and 1mM DTT.		

Introduction

Morf4 Family Associated Protein 1-Like 1 (MRFAP1L1) is a part of the MORF4 family-associated protein family.

Stability

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Background

MRFAP1L1 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 150 amino acids (1-127 a.a) and having a molecular mass of 17.2kDa. MRFAP1L1 is fused to a 23 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
